


Section 3 - How to adapt sequences and optimize images



3.1 BFFE Cine Setup



3.2 Flow quantification



3.3 Tips & tricks and Challenges



3.4 Common pitfalls and artifacts



Tips & tricks and Challenges

1. Tips & Tricks to deal with poor breath-holders(BH) and to reduce the scan time per BH with Cine sequences

Aim: follow the steps below and optimize the parameter settings.

- **Reduce number of slices acquired per breath-hold** (=increase overall scan time)
- **Increase voxel size** (=Decrease spatial resolution)
- **Use parallel imaging/acceleration technique** (Sense/CS-Sense/AI)
- On *Siemens Scanner*:
 - The **Concatenation** parameters determine the number of slices to be measured per BH. Small number of concatenation divide the measurement into only a few but long BH. A large number of concatenation divides the measurement into many short BH.

$$\frac{\text{Total number of slice}}{\text{Concatenation}} = \text{number of slice per BH}$$

- **Increase Segments** (*Siemens & GE scanners*).



Tips & tricks and Challenges

2. Tips & Tricks for CINE Free Breathing (FB)

Aim: follow the steps below and optimize the parameter settings.

➤ **Can Fetal CMR be difficult?**

- Small physique
- High heart rate
- Uncoperative with breath hold
- Fast and shallow breathing pattern

➤ ***Need manual optimization for parameter settings***

- FOV smaller
- Decrease voxel size
- Increase NSA/Average/NEX (between 2 and 5).
An increased number of signal average could help with maternal movement (respiratory)



Tips & tricks and Challenges

3. Challenges of Cine sequences

Aim: consider how important the parameters below are.

- **TR, TE and Flip Angle (FA) consideration:**
 - Keep repetition time (TR) and echo time (TE) as short as possible to maximize the image quality;
 - Avoid TR ▶ 4ms;
 - Keep Flip Angle:
 - 1.5T (70-80 is preferred)
 - 3T (60 is preferred)

- **How to further decrease the TR?**
 - TR decreases with lower resolution → increase voxel size (low matrix size & large FOV)
 - Generally, keep the frequency voxel size higher than phase voxel size






➤ **SAR**

- High SAR level are easily reached with large Flip Angle and short TR (3T)
- Reduce FA to reduce TR

➤ **Shimming box**

- Correct positioning of the shim box is essential for a good image quality and to avoid black bands artifacts & ghost artifacts from large vessels (also refer to Artifacts: How to avoid them). Therefore, always ensure to move the shim box to the region of interest



Note: Shim box volume covering the entire heart